

Dear Mr./Mrs./Ms.,


Your Doctor has prescribed you medication for your respiratory problem which is to be inhaled using the Diskus[®]. The Diskus[®] is extremely easy to use. Just open, click it, inhale and close. On the right-hand side of the page you will find detailed instructions for use.

Each Diskus[®] contains exactly 60 inhalations. Each individual one is separately packed and therefore protected against moisture. You are able to see how many single doses remain by looking in the small indicator window. The last 5 doses will appear red. The Diskus[®] is used for long-term treatment. Please stick to the recommended guidelines.

On top of that, a fast-acting preparation for the widening of the bronchial tubes is necessary for the treatment of an attack of acute shortness of breath.

1. Open

Push away thumb grip revealing mouthpiece and lever.



2. Click

Push away lever until you hear a click.



3. Inhale

Hold level. Take one deep breath in and out (not into the Diskus[®]). Place your lips firmly around the mouthpiece and breathe deeply and evenly through the Diskus[®]. Hold your breath for a short while.



4. Close

Push thumb grip back into position.

